COVID-19 & Influenza Specimen Collection and Shipment

North Carolina State Laboratory of Public Health

All COVID-19 & Influenza specimens submitted to the NCSLPH must meet the testing criteria on our website and be collected safely. The following guidance applies only to the NCSLPH.

◆ Failure to follow submission guidelines may result in specimen rejection.

Specimen Collection

The following are acceptable specimen types:
 Nasopharyngeal (NP) swab
 Oropharyngeal (OP) swab
 Nasal swab

Use a synthetic tipped, sterile swab of appropriate size with a plastic or metal shaft. Do not use a calcium alginate swab or one with a wood shaft.

- The following are acceptable media types: Viral Transport Media (VTM)
 Universal Transport Media (UTM)
- ◆ Label the vial containing the swab completely: Swab Type Patient Name ____ Or any 2 unique identifiers Date of Birth ____ Or any 2 unique identifiers
- Store properly until shipment:
 Refrigerate at 2-8° C for up to 72 hours
 Freeze ≤ -70° C for longer storage
- Completely fill out the COVID-19 submittal form online in our <u>ETOR portal</u>.

Specimen Shipment

- Specimen collection supplies and packaging and shipping supplies can be ordered online from NCSLPH Online Supply Ordering System
- All specimens shipped via commercial courier must meet <u>Category B, UN3373</u>
 requirements



- ◆ Keep cold. Ship refrigerated specimens on frozen cold packs and any ≤ -70°C frozen specimens on dry ice (the latter following IATA regulations). If shipping using a medical courier, no cold packs are required but specimens must be refrigerated.
- Label the package completely:

Attention: Virology/Serology Unit North Carolina State Laboratory of Public Health 4312 District Drive Raleigh, NC 27607-5490

 Specimen deliveries are only received at the State Lab loading dock from:

> 8 a.m. - 5 p.m., Monday - Friday 8 a.m. - noon, Saturday

Result Reporting

Results are posted electronically to our <u>CELR</u> and <u>ETOR</u> systems, tied to the submitter's EIN number. Visit our <u>website</u> for account setup and tutorials.

